


So... you broke some glass...

Safety Moment

9/19/14



Handling glassware

- Always check glassware for cracks or chips
 - Don't use broken glassware
- Take care handling glass from the oven
 - Hot glass looks like cold glass
- Don't rapidly heat or cool glassware



Oh no!

**you broke some
glassware!**

**WHAT SHOULD
YOU DO NOW?!**



First, check for
any cuts on
yourself



Next, get a broom and dust pan to clean up the glass



**Finally, place
the glass in
the broken
glass disposal**



A few more things...

- Once glass disposal is full, seal with tape and place outside of the lab.
 - It will get picked up and you will get a new one shortly
- Be sure to secure all your glassware properly while using it.
 - Clamps, ring stands, cork stands
- Also be sure not to clutter up any area with too much glassware